
FORA TNG

Blood Glucose Monitoring System

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Dear Owner of FORA TN'G Blood Glucose Monitoring System:

The system consists of three main products: the meter, test strips and control solutions. These products have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Only use **FORA TN'G** test strips and **FORA** control solutions with the **FORA TN'G** Blood Glucose Monitoring System.

TABLE OF CONTENTS

IMPORTANT SAFETY INSTRUCTIONS	04
Read Before Use	04
Important Information	06
Intended Use	07
Test Principle	07
Contents of System	08
Appearance and Key Functions of The Meter	09
Meter Display	10
Test Strip	11
Setting the Meter and Deleting the Memory	12
CHECKING THE SYSTEM WITH FORA CONTROL SOLUTIONS	15
How to Perform A Control Solution Test	15
TESTING YOUR BLOOD	18
The Four Measuring Modes	18
How to Perform A Blood Glucose Test	19
Preparing the Puncture Site	20
Performing A Blood Glucose Test	21
VIEWING THE METER MEMORY	23
Viewing Results	23
Viewing Day Average Results	24
DOWNLOADING RESULTS	26
CLEANING AND DISINFECTING YOUR METER	27
Caring for Your Meter	27
Disinfecting Procedures	28
Cleaning and Disinfecting Your Lancing Device	29
BATTERY RECHARGE	31
DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE	33
Troubleshooting	35
DETAILED INFORMATION	36
Reference Values	36
Performance Characteristics	37
SYMBOL INFORMATION	38
SPECIFICATIONS	40

IMPORTANT SAFETY INSTRUCTIONS

Read Before Use

- The meter and lancing device are for single patient use only. Do not share them with anyone including your family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
5. Before using this device to test blood glucose, read all instructions thoroughly and practice the test procedure. Carry out all the quality control checks as directed.
6. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
7. Use of this device in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
8. Do **NOT** use this device in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.
9. Do **NOT** under any circumstances use the device on newborns or infants.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to inaccurate results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients who are in shock. Readings which are lower than actual values may occur for individuals in a hyperglycemic-hyperosmolar state, with or without ketosis. Please consult your healthcare professional before use.

Intended Use

FORA TN'G Blood Glucose Monitoring System is intended for self testing outside the body (*in vitro* diagnostic use) by people with diabetes at home as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. This system is intended for single-patient use (lay-users at home) and should not be shared.

It should not be used for the diagnosis of or screening for diabetes mellitus, or testing on newborns.

This system provides you with plasma equivalent results.


Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Contents of System

Your new Blood Glucose Monitoring System includes:

- 1 Meter
- 2 1 FORA Lancing Device
- 3 10 FORA Lancets
- 4 Owner's manual
- 5 Storage Case
- 6 Quick Start User Guide / Warranty Card
- 7 Lancing Device Insert
- 8 Daily log book


Test strips, control solutions, or sterile lancets are **REQUIRED** for blood glucose testing, but may not be included in the kit (please check the contents on your product box). They can be purchased separately. Make sure you have those items needed for a blood glucose test beforehand. Please contact the retail pharmacies (the place of purchase) for availability.

NOTE:

Check your system to ensure that it is unopened prior to use and that it contains all parts listed above. If anything is missing or damaged, please return your system to the place of purchase.

Appearance and Key Functions of The Meter


- 1 Display Screen**
- 2 Test Slot**
- 3 Main Button**
- 4 Charging Indicator/Bluetooth Indicator**
- 5 SET Button**
- 6 Battery Cover**
- 7 USB Port (Only for battery recharge)**

NOTE:

The meter will turn off automatically after 180 seconds of inactivity, or you can press and hold the main button for 3 seconds to turn off the meter.

Meter Display

1 Error Message Symbol

2 Measurement Unit

3 Testing Result Display

4 Measuring Mode

Gen, AC, PC, QC.

5 Memory Symbol

Indicates a test result stored in memory.

6 Blood Drop Symbol

Flashes when meter is ready for sample application.


7 Time Display

8 Alarm Symbol

9 Low Battery Symbol

Test Strip

Your system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole on the test strip and is automatically drawn into the reaction cell where the reaction takes place.

Attention:

The front side of test strip should face up when inserting test strip.

Test results might be wrong if the contact bar is not fully inserted into the test slot.

NOTE:

FORA TN'G Blood Glucose Monitoring System should only be used with **FORA TN'G** test strips. Using other test strips with this meter can produce inaccurate results.

Setting the Meter and Deleting the Memory

Start with the meter off (no test strip inserted). Press **SET**.

1. Setting the date

With the year flashing, press Main button until the correct year appears. Press SET.

With the month flashing, press Main button until the correct month appears. Press SET.

With the day flashing, press Main button until the correct day appears. Press SET.

2. Setting the time format

Press Main button to select the desired time format — 12h or 24h. Press SET.

After you have finished selecting the 12h or 24h clock, press SET and then the hour will flash on the screen.

3. Set the hour and minute

With the hour flashing, press Main button until the correct hour appears. Press SET. With the minute flashing, press Main button until the correct minute appears. Press SET.

4. Deleting the memory

With “dEL” and a flashing “**M**” on the display, press Main button and select “no” to keep the results in memory then press SET to skip.

To delete all the results, press Main button and “yes” and “**M**” are displayed on the meter, press SET to delete the memory.

If you do not want to delete the results, do not press the Main button but press the SET button again to turn off the meter.

5. Set the reminder alarm

You may set up any or all of the reminder alarms (1-4). The meter displays “On” or “OFF” and “”, press Main button to turn on or turn off to set the first reminder alarm.

Press Main button to select “On”, then press SET to set the hour. When the hour is flashing, press Main button to add an hour. Press SET to confirm and go to minutes, press Main button to add one minute. Hold Main button longer to add

faster. Press SET to confirm and go to the next alarm setting.

If you do not want to set an alarm, press SET to skip this step.

If you want to turn off an alarm, find the alarm number by pressing SET in the setting mode, press Main button to change from “ON” to “OFF”.

6. Enter the Bluetooth pairing

With “no” flashing on the display, press SET to skip the pairing.

If you wish to enter the pairing mode, press Main button once and the meter will display “yes”. Then press SET to confirm and the meter will show “CLr”. The meter will then automatically turn off and the Bluetooth indicator will start blinking, this means the meter is in the process of pairing with a Bluetooth receiver.

Another way to enable Bluetooth pairing is to press and hold Main button in memory mode (refer to chapter VIEWING THE METER MEMORY) until the meter turns off. While the Bluetooth indicator blinks, you can pair the meter with a Bluetooth receiver.

NOTE:

This step is recommended when the user needs to pair the meter to a Bluetooth receiver for the first time, or when the user needs to pair the meter to another new Bluetooth receiver.

CHECKING THE SYSTEM WITH FORA CONTROL SOLUTIONS

FORA control solutions contain a known amount of glucose that reacts with test strips. By comparing the result of your control solution test with the expected range printed on the test strip vial label, you can check if the meter and test strips are working properly, and if you're performing the test correctly. It is very important that you perform this simple check routinely to make sure you get accurate results.

Do a control solution test when:

- you first receive the meter,
- at least once a week to routinely check the meter and test strips,
- you begin using a new vial of test strips,
- you suspect the meter or test strips are not working properly,
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
- practicing the testing process, or
- you dropped or think you have damaged the meter.


How to Perform A Control Solution Test

FORA control solutions contain a known amount of glucose that reacts with test strips and is used to ensure that the meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have all items prior to testing your blood glucose. Please contact the Customer Care Line at 1-888-307-8188 or the place of purchase for availability.

1. Insert test strip


Insert a test strip into the test slot with the contact bars end first and facing up. (Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.) The meter turns on automatically and displays the following in sequence:

When the “” appears on the display, press the Main button and “QC” will appear on the display. When the “QC” sign is displayed, the meter will not store your test result in the memory. If you decide not to perform a control solution test, press the Main button again and the “QC” sign will disappear.

CAUTION:

Every time you perform a control solution test you must enter “QC” test mode so that the test result will not be stored in the meter memory.

2. Apply FORA control solution

Shake **FORA** control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter in an angled position so that the absorbent hole of the test strip touches the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply **FORA** control solution onto the strip.

3. Read and compare the results

After the meter counts to 0, the result of the control solution test will appear. Compare the result with the range printed on the test strip vial. The result should fall within this range.

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Contact customer service or place of purchase for help.

NOTE:

- The control solution range printed on the test strip vial is for **FORA** control solution use only. It is not a recommended range for your blood glucose level.
- See the Meter Maintenance section for important information about your FORA control solutions.

TESTING YOUR BLOOD

Before testing, be sure to read this section and the test strip insert found in the test strip box carefully. Make sure you have all the items needed for the test.

The Four Measuring Modes

The meter provides you with four modes for measuring: General, AC, PC and QC.

MODES	USE WHEN
General tests	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution

You can switch between the four modes by:

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a “”.

2. Press Main button to switch between General, AC, PC and QC mode.

How to Perform A Blood Glucose Test

Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. The lancing device must not be used for more than one person.

To reduce the chance of infection:

- **Wash and dry your hands before starting.**
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Wash and dry your hands thoroughly after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
- If the meter is being operated by a second person who is providing

testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.

Preparing the Puncture Site

Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site on fingertips.
- Rub the puncture site for about 20 seconds before penetration.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Use a clear cap while setting up the lancing device.

• **Fingertip testing**

Press the tip of the lancing device firmly against the lower side of your fingertip. Press the release button to puncture your finger, a click indicates that the puncture is complete.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Performing A Blood Glucose Test


1. Insert the test strip to turn on the meter

Wait for the meter to display “”.


2. Select the appropriate measuring mode by pressing Main button

For selecting the measurement mode, please refer to “The FOUR MEASURING MODES”.

3. Apply blood sample

Obtain a drop of blood of at least 0.5 μL with a lancing device. When the “” symbol appears on the display, touch the drop of blood with the absorbent hole of the test strip. Do not rub the test strip with your finger or try to apply a smeared sample.

Gently apply the drop of blood to the absorbent hole of the test strip at a tilted angle. Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a beep sound.

If the confirmation window is not filled completely before the meter begins to count down, do not add more blood to the test strip. Discard the test strip and start again. If you have trouble filling the test strip, please call your local customer service number for assistance.

NOTE:

If you do not apply a blood sample within 3 minutes, the meter will automatically turn off. You must remove the test strip and insert it back into the meter to restart the test.

**4. Get result**

The result of your blood glucose test will appear after the meter counts down to 0. The blood glucose result will be stored in the memory automatically.

5. Eject the used test strip

Remove the test strip by hand, throw the strip into a disposal container. The meter will turn itself off automatically after the test strip is removed.

Always follow the instructions in the lancing device insert when removing the lancet.

WARNING:

- Please do not change your treatment based on the result without first consulting your health care professional.
- Turn the meter off by removing the test strip. Discard the used test strip and lancet carefully according to your local regulations.
- Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination.

VIEWING THE METER MEMORY

Your Meter stores the 450 most recent blood glucose test results with date and time in its memory. You can review the test results with these easy steps.

Viewing Results

1. Enter the memory mode

With the meter turned off, press the Main button. The first test result will appear, indicating that you are in the memory mode.


When using the meter for the first time or if the results have been deleted, “- - -” will appear, indicating that there are no test results in the memory.

2. Recall test results

After the last test result, the most recent test result appears with the date and time. Press the Main button to review your last 450 test results in order. When the memory is full, the oldest result will be deleted when the newest is added.

3. Exit the memory mode

Press the Main button for three seconds to turn off the meter or after the last result, the display will show “End” and it will turn off automatically.

NOTE:

The results of control solution tests are **NOT** stored in the memory. Stored results are blood glucose results only.

Viewing Day Average Results

1. Press and release Main button

Keep pressing Main button for 3 seconds until the flashing “AVG” appears. Release Main button and then your 7-day average result measured in general mode will appear on the display.

2. Press Main button to review

14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.

3. Exit the meter memory

Keep pressing the Main button and the meter will turn off after displaying the last test result.

Note:

- Any time you wish to exit the memory, keep pressing Main button for 5 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are NOT included in the day average.
- If using the meter for the first time, “--” displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.

DOWNLOADING RESULTS

Data Transmission via Bluetooth

You can use your device with iOS or Android system to download data from your meter via Bluetooth. Follow the steps below to transmit data from your meter. Please contact local customer service or place of purchase for assistance.

Please note that you must complete the pairing between meter and Bluetooth receiver before transmitting data.

- 1. Make sure your meter is already paired with your device with iOS or Android system by following instructions of the setting mode in Step 6.**
- 2. Install the software on your device with iOS or Android system.**
- 3. Data transmission**

Turn on the Bluetooth on your device with iOS or Android system and keep it on, and your meter will be able to connect to the device after measurement. You can view the data on your device.

Bluetooth indicator on the meter

BLUETOOTH INDICATOR	STATUS
Flashing Blue	The Bluetooth function is enabled and waiting for connection.
Solid Blue	The Bluetooth connection is established.

NOTE

- While the meter is in transmission mode, it will be unable to perform a test.
- Make sure your iOS (6 or above) or Android system (4.3 or above) device has its Bluetooth enabled and your meter is within receiving range before you transmit data from your meter.
- A compatibility issue between your mobile device and the meter might occur due to the way Bluetooth functionality is implemented by the various mobile device manufacturers.
- The USB port of the meter is only for battery recharge. It does not have the function of data transmission.

CLEANING AND DISINFECTING YOUR METER

Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

Why the cleaning and disinfection should be performed

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

How to clean and disinfect the meter

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. Do **NOT** use organic solvents to clean the meter.

We recommend for meter cleaning and disinfection you should use the disinfecting wipe from below, which the active ingredients have been tested to be effective against Hepatitis B Virus (HBV) for the meter.

► Clorox Healthcare™ Professional Disinfecting Bleach Wipes
(DIN no: 02401983)

The wipe has been shown to be safe for use with the meter.

Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

2. Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep the meter surface wet with disinfection solution for a minimum of 1 minutes for Clorox Healthcare™ wipes. Follow the instructions on the package label of disinfecting wipe.

3. Allow the meter surface to dry completely.

4. Discard the used wipes and never reuse them.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at 1-888-307-8188 for assistance.

This device has been validated to withstand 260 cleaning and disinfection cycles, using the recommended disinfecting wipe/towelette, to simulate one cleaning and disinfection cycle per week over 5 years of use.

The meter should be replaced after the validated number of cleaning and disinfection cycles or the warranty period, whichever comes first.

Stop using the meter if you see any signs of deterioration. For example:

- meter can not be turned on,
- LCD display cracks or becomes cloudy,
- buttons no longer function,
- meter outer casing cracks,
- data can not be transmitted to PC,
- color or paint/printing on housing is abnormal, or
- scratches or abrasions on meter are higher than acceptable.

Please contact the customer service for a replacement meter if any of the signs of deterioration are noticed.

Cleaning and Disinfecting Your Lancing Device

The lancing device is reusable. Please disinfect your lancing device regularly.

- You should disinfect the lancing device at least once a week to prevent infection. Or clean the lancing device when you see any dirt on it.
- To clean the lancing device, use a disinfecting wipe to clean the lancing device. Please refer to disinfection procedure for more details.
- Do NOT place the device in a dishwasher or use detergents.
- We recommend for lancing device cleaning and disinfection you should use the disinfecting wipe from below:
- Clorox Healthcare™ Professional Disinfecting Bleach Wipes (DIN no: 02401983)

NOTE :

- Do **NOT** clean and disinfect the meter while performing tests.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person. and the user must wear gloves during cleaning and disinfection.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot.

- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

BATTERY RECHARGE

Your meter comes with a built-in rechargeable Li-polymer battery. If the low battery symbol “” appears on the screen indicating that the battery is low and it is time to recharge the battery. We advise you to fully recharge every 9 months when not in-used to ensure best condition.

1. Connect the USB cable to the data port of the meter.
2. Connect the other end of cable to a free USB port on your personal computer.
3. “USB” will appear on the display and the charging indicator will light up in red when battery is recharging. After the battery is fully charged, the red light will turn green.
4. Remove the USB cable and the meter will automatically switch off.

NOTE :

- Meters may not be rechargeable if recharging is not performed over 9 months.
- It should take approximately 2 hours to be fully recharged. Recharging the battery does not affect the test results stored in the memory.
- The USB cable is for battery recharging and data transmission with a computer. If you wish to recharge through a wall outlet, please obtain a USB adapter.

Meter Storage:

Storage conditions: -20 °C to 60 °C, below 95% relative humidity.

- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage conditions: 2 °C to 32 °C, below 85% relative humidity. Do NOT freeze.

- Store your test strips in their original vial only. Do not transfer to another container.
 - Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
 - After removing a test strip from the vial, immediately close the vial cap tightly.
 - Touch the test strip with clean and dry hands.
 - Use each test strip immediately after removing it from the vial.
 - Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months.
 - Do not use test strips beyond the expiry date. This may cause inaccurate results.
 - Do not bend, cut, or alter a test strip in any way.
 - Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.
- For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only **FORA** control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 20 °C to 25 °C. Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2 °C to 30 °C. Do NOT freeze.

DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE

The following is a summary of display messages. If your meter displays an error message, please follow the actions for the error message as described in the table below. If the problem persists, please contact your local customer service agent for help.

MESSAGE	WHAT IT MEANS	ACTION
	Appears when glucose < 1.1 mmol/L	This indicates hypoglycemia (low blood glucose.) Please seek immediate medical assistance.
	Appears when glucose > 33.3 mmol/L	This indicates severe hyperglycemia (high blood glucose). Please seek immediate medical assistance.
	The “  ” symbol appears with this message. This means that the battery does not have enough power remaining for a test.	Replace the battery immediately.
	A used test strip is inserted, or there is an electrical problem.	Repeat test with a new test strip. If the error message re-appears, please contact your local customer service agent for help.
	You may have removed the strip after applying blood to the absorbent hole.	Keep the test strip in the meter after it has absorbed the blood until the meter displays the test result.

MESSAGE	WHAT IT MEANS	ACTION
	<p>The temperature of the environment, meter, or test strip is outside the system's operating range. You cannot perform a test until the meter and test strip are within the operating range of 10 °C to 40 °C</p>	<p>Repeat the test after the meter and test strip are within the operating temperatures.</p>
	<p>Problem with the meter.</p>	<p>Review the instructions and try again with a new test strip. If the problem persists, please contact your local customer service agent for help.</p>

For questionable results, please see test strip package insert for important information.

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer service.

2. If the test does not start counting down after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer service.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake FORA control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of FORA control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature (20 °C to 25 °C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer service.

DETAILED INFORMATION

Reference Values

Blood Glucose

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that maintaining **blood glucose levels close to normal** can reduce the risk of diabetes complications by up to 60%*¹. The results provided by this system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 5.6 mmol/L
2 hours after meals	< 7.8 mmol/L

Source: American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.

Please work with your doctor to determine a target range that works best for you.

References:

*1: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

PERFORMANCE CHARACTERISTICS

Accuracy

Within ± 0.83 mmol/L at glucose concentration < 5.5 mmol/L and within $\pm 15\%$ at glucose concentration ≥ 5.5 mmol/L. Studies were conducted using capillary and the results are as below:

N= 160
$y = 0.9713x + 3.8502$
$R^2=0.9795$
Range=2.26-30.0 mmol/L

The tables show differences in glucose values between this system and YSI-2300.

Table 1 represents samples for glucose results lower than ($<$) 5.5mmol/L

Table 2 represents samples for glucose results greater than (\geq) 5.5mmol/L.

Table 1 Glucose results lower than ($<$) 5.5 mmol/L

Difference range in values between the laboratory value and FORA TN'G's value:	Within 0.28 mmol/L	Within 0.56 mmol/L	Within 0.83 mmol/L*
The percent (and number) of samples for which the difference between FORA TN'G's and laboratory value were the difference range shown in the row:	181/300 (60.3%)	273/300 (91.0%)	298/300 (99.3%)

Table 2 Glucose results greater than and equal to (\geq) 5.5 mmol/L

Table 2 Glucose results greater than and equal to (\geq) 5.5 mmol/L	Within 5%	Within 10%	Within 15%
The percent (and number) of samples for which the difference between FORA TN'G's and laboratory value were the difference range shown in the row:	317/660 (48.0%)	557/660 (84.4%)	641/660 (97.1%)

*Acceptance criteria in ISO 15197:2013(E) are that 95% of all differences in glucose values (i.e., YSI-2300 glucose value minus FORA TN'G's glucose value) should be within ± 0.83 mmol/L of an YSI-2300 laboratory method for glucose values less than 5.5 mmol/L, and within 15% for glucose values greater than 5.5 mmol/L.

NOTE:

When **FORA TN'G** glucose meter results are compared to the laboratory results, difference values below 5.5 mmol/L are expressed in mmol/L, while those above 5.5 mmol/L are compared in percent.

Additional information of Accuracy for users

Diabetes experts have suggested that glucose meters should agree within 0.83 mmol/L of an YSI-2300 laboratory method when the glucose concentration is lower than 5.5 mmol/L and within 15% of a laboratory method when the glucose concentration is 5.5 mmol/L or higher. Table 3 displays how often the **FORA TN'G** meter achieves this goal. Table 3 is based on a study done on 160 patients to see how well the **FORA TN'G** compared to laboratory results.

Table 3

For glucose results lower than 5.5 mmol/L, the percent (and number) of meter results that match the laboratory method within 0.83 mmol/L:	50/50 (100%)
For glucose results at 5.5 mmol/L or higher, the percent (and number) of meter results that match the laboratory method within 15%:	107/110 (97.3%)

NOTE:

When meter results are compared to the laboratory results, results below 5.5mmol/L are compared in mmol/L.

Precision

CVs (%) of intermediate precision and repeatability were less than 5%. The product has been tested to meet the requirements of ISO 15197.

SYMBOL INFORMATION

The following symbols are used throughout the product labeling for **FORA TN'G** blood glucose monitoring system:

Symbol	Referent	Symbol	Referent
	<i>In vitro</i> diagnostic medical device		Use by
	Do not reuse		Manufacturer
	Consult instructions for use		Batch code
	Temperature limitation		Serial number
	Humidity Limitations		Caution, consult accompanying documents

SPECIFICATIONS

Identifier No.: TN'G

Dimension: 93.0 (L) x 26.0 (W) x 15.5 (H) mm

Weight: 31 g

Power Source: Li-Polymer battery

Display: LCD

Memory: 450 measurement results with date and time

External Output: Bluetooth

Automatic detection of electrode insertion

Automatic reaction time count-down

Auto switch-off after 3 minutes without action

Temperature warning

Operating Condition:

10 °C to 40 °C, below 85% R.H. (non-condensing)

Storage/Transportation Conditions:

-20 °C to 60 °C, below 95% R.H.

Measurement Units: fixed mmol/L

Measurement Range: 1.1 to 33.3 mmol/L

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, EN 61326-2-6.

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ForCare, Inc.