# FORA TN'G

## Warnings

- ▶ For *in vitro* diagnostic use (for use outside of the body only).
- ▶ For single patient use only.
- ▶ Healthcare professionals and other users testing multiple patients with this system should handle everything that comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- ▶ Please read this sheet and your Blood Glucose Monitoring System Owner's Manual before you use this test strip. **Use only FORA TN'G Test Strips with FORA TN'G Blood Glucose Meter to obtain accurate results, and be covered by the manufacturer's warranty.**
- ▶ Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- ▶ Low results may be inaccurate when testing on patients in hyperglycemic-hyposmolar state, with or without ketosis.
- ▶ Please do not use the Glucose Monitoring System on critically ill patients.
- ▶ Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

## Intended Use

FORA TN'G Blood Glucose Test Strips, when used together with FORA TN'G Blood Glucose Monitoring System, allow your blood glucose levels to be measured by yourself at home. It uses fresh whole blood samples from the finger. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

## Limitations


- ▶ Red blood cell count (Hematocrit): Your hematocrit level must be between 20% and 60%. Please ask your healthcare professional if you do not know your hematocrit level.
- ▶ Neonatal Use: **This test strip must not be used for the testing of newborns.**
- ▶ Metabolites: Dopamine, L-Dopa, methyl dopa, uric acid, ascorbic acid (vitamin C), and acetaminophen at normal blood concentration do not significantly affect blood glucose test results.
- ▶ Xylose: Do not test blood glucose during or soon after a xylose absorption test. Xylose in the blood can give falsely elevated results.
- ▶ Lipemic Effects: Blood triglycerides up to 166.7 mmol/L (3000 mg/dL) do not affect the results significantly, but may affect results at higher levels.
- ▶ Glutathione reduced and pralidoxime iodide: Do not test blood glucose during or soon after a glutathione reduced or pralidoxime iodide treatment. The two compounds within the therapeutic or physiologic concentration range may affect the glucose results.
- ▶ Bilirubin, tolazamide, galactose, mannitol, and mannose may affect the glucose results at high levels.
- ▶ Altitude Effects: Altitudes up to 3,275m (10,742 feet) do not affect test results.

## Storage and Handling

**IMPORTANT:** Do not use the test strips if they have expired.

- ▶ Test strips expire 3 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- ▶ Store the test strips in a cool, dry place between 35.6°F and 89.6°F (2°C and 32°C) and 10% to 85% relative humidity.
- ▶ Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- ▶ Store the test strips in their original vial **ONLY**. Do not transfer them to a new vial or any other containers.
- ▶ Do not touch the test strips with wet hands.
- ▶ Use each test strip immediately after taking it out of the vial. Close the vial immediately after taking out a strip.
- ▶ Keep the vial closed at all times.
- ▶ Do not bend, cut, or alter the test strip.

## Strip Appearance

### 1. Absorbent Hole

Apply a drop of blood here. The blood will be automatically absorbed.

### 2. Confirmation Window


This is where you confirm if enough blood has been drawn into the absorbent hole of the strip.

### 3. Test Strip Handle

Hold this part to insert the test strip into the slot.

### 4. Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.


Front

**ATTENTION:** The front side of test strip should face up when inserting test strip. Test results might be wrong if the contact bar is

not **fully** inserted into the test slot.


## Testing Your Blood Glucose

**PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.**

### STEP 1

Insert the test strip fully into the slot of the meter until it will go no further. When the strip is fully inserted, the meter will do several self-checks.

### STEP 2

Collect a blood sample with the test strip. A sufficient quantity of blood is required for the test to provide accurate results. Touch the blood drop with the absorbent hole of the test strip, and wait until the confirmation window is fully covered. **Do NOT** apply a smeared blood sample. The meter will start counting down.

5.1  
mmol/L

### STEP 3

After a few seconds, the meter will display your blood glucose level. The last reading will be automatically saved in the meter. Turn it off by removing the test strip and throw away the used test strip.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations

## Reading Your Result

Your blood glucose readings deliver plasma equivalent results and are displayed in millimoles of glucose per litre of blood (mmol/L).

### Reference values

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meals	Less than 5.6 mmol/L
2 hours after a meal	Less than 7.8 mmol/L

Source: American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.

**Please consult your doctor to determine a target range that works best for you.**

### Questionable or inconsistent results

If your test results are unusual or inconsistent with how you are feeling:

- Make sure the confirmation window of the test strip is completely filled with blood.
- Check the expiration date of the test strips.
- Check the performance of your meter and test strip with the control solutions.

**PLEASE NOTE:** Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

## Quality Control Testing

**IMPORTANT:** The reference range of the control solutions may vary with each new vial of test strips. Make sure you check the range on the label of your current vial.

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

## Chemical Components

> Glucose dehydrogenase (E. coli)	8%
> Electron shuttle	55%
> Enzyme protector	8%
> Non-reactive ingredients	29%

## Performance Characteristics

### Accuracy

Within  $\pm 0.83$  mmol/L at glucose concentration < 5.5 mmol/L and within  $\pm 15\%$  at glucose concentration  $\geq 5.5$  mmol/L.

### Precision

CVs (%) of intermediate precision and repeatability are less than 5%. The product has been tested to meet the requirements of ISO 15197.

## Symbol Information

Symbol	Referent
IVD	For <i>In vitro</i> diagnostic use only
Ⓜ	Do not reuse
i	Read instructions before use
Ⓜ	Humidity limitations
Ⓜ	Temperature limitations

Symbol	Referent
Manufacturer	Manufacturer
⚠	Caution, consult accompanying documents
SN	Serial number
📅	Use by
LOT	Batch code

### ForaCare, Inc.

893 Patriot Dr., Suite D, Moorpark, CA 93021 USA

### ForaCare Customer Care:

**1-888-307-8188** (7:00am-6:00pm PST, Mon.-Fri.)

For assistance outside of these hours, please contact your healthcare professional. [www.foracare.com](http://www.foracare.com)

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